## Viscosity Cups - Ford, DIN, ISO

Operating Manual

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# Table of Contents

English .....	4
Deutsch.....	8

# 1 Description

Viscosity cups have been produced by BYK Gardner for over forty years. The product range comprises cups for following standards:

- DIN 53211\* (was withdrawn in October 1996)
- DIN EN ISO 2431, ASTM D 5125
- ASTM D 333, D 365, D 1200

For many applications it is not necessary to know the absolute viscosity of a paint or varnish. A parameter permitting a relative classification and estimation is considered sufficient. The efflux time, measured in seconds, is such a reliable parameter. It is determined using flow cups of various designs. These cups hold a defined volume of liquid which flows through an orifice.

The reproducibility of such measurements depends on:

- The accuracy of the size of the cup,
- a constant temperature during measurement,
- the Newtonian flow of the liquid.

BYK-Gardner guarantees the accuracy of the size of its flow cups by using NC machines for the manufacturing process, surface refinement, quality control of raw materials through the manufacturing process until shipment, and by examining every single piece using test oil in a climatic chamber. Sufficiently long tempering of the cup and test sample produces a constant test temperature. For measurement at temperatures differing from room temperature, the cup is put into a temperature control jacket, which contains a liquid that is heated up or cooled down by means of a thermostat. The liquid under test flows over a thermometer, thus allowing the actual test temperature to be determined.

For efflux time measurement it is important to note that a cup can only produce satisfying and reproducible results when used within its defined viscosity range.

## **NOTICE**

Flow cups cannot measure high viscosity liquids, pastes etc.

## 2 Operation

Viscosity cups are used to determine the viscosity of a sample by measuring the efflux time of the test material.

### Measurement Procedure

Perform following steps to determine the efflux time:

1. Temper the cup and test liquid to e.g. 20 °C ± 0.5 °C (68 °F ± 1.0 °F).
2. Put flow cup in a precisely horizontal position into the ring stand or temperature control jacket.
3. Close orifice.
4. Pour in test liquid.
5. Draw clean glass plate over the rim of the cup, thus removing superfluous liquid into overflow reservoir and closing the cup.
6. Open orifice.
7. Remove glass plate horizontally and start stopwatch.
8. Stop the watch with the first break in the efflux stream.
9. Repeat three times each with a new sample of the same material.

### Conversion Table

The efflux time in seconds corresponds to the kinematic viscosity in centistokes (cP). For this a BYK-Gardner provides a "Viscosity Cup Conversion Chart".

7201 0172	0175 0173	0176 0174	0115 0115	0213 0213	0214 0214	0215 0215	0216 0216	Part No. with cert
Ford Cups			DIN Cup 4mm	ISO Cups				
2	3	4	4mm	3mm	4mm	5mm	6mm	
Drain Time in sec			Drain Time in sec	Drain Time in Sec				Viscosity In Centistokes (mm <sup>2</sup> /s)
								2.0
								2.5
								3.0
								3.5
								4.0
								4.5
								5.0
								5.5
								6.0
								6.5
								7.0
								7.5
								8.0
								8.5
								9.0
								9.5
								10.0
								10.5
								11.0
								11.5
								12.0
								12.5
								13.0
								13.5

Select the viscosity cup used, search for the measured efflux time in seconds and depict the viscosity in centistokes (cP) on the right side.

#### NOTICE

The conversion chart is provided as separate document. The document number is "199 023 029". This chart can be downloaded along with the operation manual, e.g. under: <https://www.byk-instruments.com/p/0115>. The chart covers all supported cups.

## 3 Cup Details

### Ford Cup in accordance with ASTM D 1200, D 333, D 365

Each Ford cup produced has been in conformance with viscosity cup flow characteristics as initially listed in ASTM Bulletin No. 169 of October 1950. The Ford cups produced by BYK-Gardner comply with the requirements listed in ASTM D 1200.

The ASTM D 1200, D 333, D 365 standard prescribes 3 cups with different orifices:

- Ford Cup No. 2 > 2,53 mm vis. at 25 °C = 25 - 120 mm/s
- Ford Cup No. 3 > 3,40 mm vis. at 25 °C = 40 - 220 mm/s
- Ford Cup No. 4 > 4,12 mm vis. at 25 °C = 70 - 370 mm/s

For measuring procedure see description of DIN cup. The efflux times measured should not fall below 20 s or exceed 100 s.

### Flow Cup in Accordance with DIN 53211 \*

This cup holds 100 ml ± 1 ml and has an integrated orifice with an orifice diameter of 4 mm ± 0.02 mm. The cup is made of anodized aluminum while the orifice is made of stainless steel, with the interior polished.

Generally, this DIN cup is used for liquids of viscosities ranging from approx. 100 - 500 cSt (mm<sup>2</sup>/s), which corresponds to an efflux time of 20 - 110 s.

### Flow Cup in Accordance with DIN EN ISO 2431, ASTM D 5125

The international standard describes a process to determine the efflux time using a 4 mm cup with dimensions differing from the DIN 53211\* cup. A longer orifice, less tapered body and slightly different inner dimensions provide efflux times different from those of DIN 53211\* cups, thus extending the measurement range and making the DIN EN ISO cup a useful supplement to the DIN cup. For measuring procedure see description of DIN cup. DIN EN ISO cups are primarily recommended for international use.

In addition to the 4 mm orifice, the DIN EN ISO standard describes 3-, 5- and 6-mm orifices. The most favorable measuring range for each cup is:

- DIN EN ISO Cup 3 mm > viscosity 10 - 40 mm/s = 30 - 100 s
- DIN EN ISO Cup 4 mm > viscosity 25 - 130 mm/s = 25 - 100 s
- DIN EN ISO Cup 5 mm > viscosity 70 - 360 mm/s = 25 - 100 s
- DIN EN ISO Cup 6 mm > viscosity 130 - 700 mm/s = 25 - 100 s

\* DIN 53211 was withdrawn in October 1996.

## 4 Maintenance

### Recalibration

The cups are milled from solid aluminum bar stock and are fitted with a stainless-steel orifice. Each cup is tested with a viscous fluid referenced to National Bureau of Standards certified oils. The orifice is honed to ensure that the assembly complies to the applicable "time viscosity" flow formula. If the orifice of the cup becomes damaged or worn and the original calibration is questioned, the cup should be returned to BYK-Gardner, to be checked, fitted with a new orifice (if necessary) and recalibrated.

### Care and Maintenance

Clean the cup by use of a suitable solvent and a cleaning swab or soft brush following each use. In cleaning the orifice, take special care to avoid leaving any film deposits on the inside wall. Under no condition should cleaning tools made of metal be brought into contact with the orifice.

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**BYK-Gardner GmbH**